

# Revolutionizing the approach to cardiorespiratory fitness

Presentation on VentriJect and the Seismofit® solution





# Cardiorespiratory fitness

- Cardiorespiratory fitness (CRF) represents the body's maximum ability to absorb oxygen and is a measure of one's physical ability and overall health. It is quantified as VO<sub>2</sub>max.
- CRF is mainly determined by the function and size of the heart.
- The importance of CRF is strongly acknowledged, and several studies are showing a clear link between CRF and excess mortality rates.

# Importance of cardiorespiratory fitness



Cardiorespiratory fitness (CRF), quantifiable as  $VO_2\text{max}$ , should be regularly examined as a clinical vital sign in general clinical practice as probably the **strongest indicator of health.**

American Heart Association (AHA)



We (again) call on both clinicians and public health professionals to adopt Cardiorespiratory fitness (CRF) as a **key health indicator.**

Editorial comment, Journal of the American College of Cardiology



# Importance of cardiorespiratory fitness

## Higher mortality

A registry study of 122,000 patients concluded that the low fitness group had a **5 times higher mortality** versus the high fitness group over 10 years.

Paraphrased from Mandsager, Kyle, et al: JAMA network open 1.6 (2018): e183605-e183605.

## Risk of cancer in men

If you have a high VO<sub>2</sub>max there is a **significant lower risk of colon and lung cancer** for men, as well as a **lower risk of dying from cancers** in the colon, lung and prostate.

Paraphrased from Ekblom-Bak et al; JAMA Network Open 2023

## Relatively small changes have significant impact

A decline in CRF of >2.0 METS was associated with a 74% increase in risk for low-fit individuals with Cardiovascular Diseases (CVD), and 69% increase for those without CVD.

Kokkinos et al.: J Am Coll Cardiol. 2023 Mar, 81 (12) 1137–1147

# The impact of training

Small increases in Cardiorespiratory Fitness (CRF) (eg 3-6 VO<sub>2</sub>max-units) are associated with **10-30% reduction in adverse cardiovascular events.**

Kodama et al: JAMA. 2009; 301:2024-2035

Exercise performed frequently over weeks or months may **improve CRF with 15-20%** in VO<sub>2</sub>max in adults.

Skinner et al.:Med Sci Sports Exerc. 2000; 32: 157-161 (HERITAGE)



## Examples of current approaches for estimating $\text{VO}_2\text{max}$




CPX / CPET  
(‘the gold standard’)



Submaximal  
ergometer test



Chester  
step test



Sport watches  
(Garmin, Polar, etc.)

# Current approaches have many flaws

- They all require physical exertion.
- They are all **time-consuming**.
- Most of them are **imprecise**.
- **Regular calibration** will be needed for some of the approaches.
- Some require **large expensive equipment** that is **difficult to transport**.
- **Highly trained staff** is required for some of the approaches.

# The solution: Seismofit<sup>®</sup>

An easy-to-use, scalable and precise solution to estimate VO<sub>2</sub> max

- ✓ Patient will be resting during test (no physical exertion)
- ✓ Takes less than 3 minutes
- ✓ No calibration
- ✓ Very limited education required



The patient goes to health check



The patient lies down with Seismofit® for 40 seconds



Seismofit® data uploads to the cloud



Data gets processed by Seismofit® algorithm




Result shows on smartphone



All done in less than 3 minutes



# Small, but mighty




CE-marked device  
(class 1 medical device)

Patented technology

Clinical validated


# Utilizing SeismoCardioGraphy (SCG)



SCG is accelerometer recordings of chest wall vibrations caused by the beating heart.



First demonstrated in 1957.

SCG research has bloomed in recent years due to MEMS accelerometers and clinical break through is underway.





Overall mean electrocardiogram and seismocardiogram signal . The circles indicate the mean location of the following physiologic events found in ultrasound images, from left: Atrial systole (AS), peak atrial inflow (PAI), mitral valve closure (MC), aortic valve opening (AO), peak systolic outflow (PSO), aortic valve closing (AC), mitral valve opening (MO), early ventricular filling (EVF). The grey areas indicate the 95% confidence intervals of the means for the physiologic events found in the ultrasound images \*.

\*Sørensen, K., Schmidt, S.E., Jensen, A.S. et al. Definition of Fiducial Points in the Normal Seismocardiogram. *Sci Rep* 8, 15455 (2018).  
 \*\*R. S. Crow et al., "Relationship between seismocardiogram and echocardiogram for events in the cardiac cycle," *American Journal of Noninvasive Cardiology*, vol. 8, pp. 39–46, Jan. 1994.  
 \*\*\*Agam, Ahmad, et al. "Correlation between diastolic seismocardiography variables and echocardiography variables." *European Heart Journal-Digital Health* 3.3 (2022): 465-472.

# SCG is a measure of cardiac function

- Systolic and diastolic events such as heart valve opening and closing are reflected in the SCG signal.
- Correlation between SCG fiducial points and hemodynamic events has been validated using echocardiography.

# Determining $VO_2\text{max}$

VentriJect uses AI and an advanced algorithm for determining  $VO_2\text{max}$ .

VentriJect's algorithm for prediction of  $VO_2\text{max}$  is based on:

Basic physical  
attributes

Gender

Height

Age

Weight

SCG  
measurements


Timings

Morphology

Amplitudes

Variability

Frequencies




# Seismofit<sup>®</sup> System

## Seismofit<sup>®</sup> Sensor



Seismofit<sup>®</sup> measures seismocardiography vibrations arising from heartbeats and transmitted to the chest. Each heart has its own signature relating to the opening and closing of valves and the relaxation and contraction of the heart chambers.

## Seismofit<sup>®</sup> App



The App for the mobile controls seismocardiography recordings and can be downloaded from App-stores. The App leads the user securely through a recording by instructive screens and images.

## Seismofit<sup>®</sup> Patch



Patches are for fixation of the sensor to the chest at the distal part of the sternum of the test person. The patch consists of hypoallergenic adhesive to both surfaces covered with release liners.

## Seismofit<sup>®</sup> Algorithm




The Algorithm on the VJ-server converts the recording into an equivalent to  $VO_2$ max score by use of features from the heart like timing, morphology, amplitude, peak to peak ratios and HRV in the recorded spectrum of the heartbeats.


# VentriJect Architecture

## Cloud based infrastructure

- Ensures a constant data flow
- Rapid updates to software and algorithms
- Pay per service




Example of the same patient measured on 3 different days





# Seismofit<sup>®</sup> has high reproducibility

Seismofit<sup>®</sup> is at least as reproducible as today's gold standard of  $VO_2max$  measurement.



Overall better performance in test-set  
(n=98)

MAPE (Mean absolute percentage error) compared to other tests



*Seismofit® performance from Schmidt, S. E., Hansen, M. T., Sørensen, K., Rømer, T., Gliemann, L., Karbing, D. S., Poulsen, M. K., Helge, J. W., & Søgaard, P. (2022). A Chest-Mounted Accelerometer for Estimation of Cardiorespiratory Fitness. Cardiovascular Digital Health Journal, 3(4)*

# Comparison to other test methods

- Seismofit® comes close to the Gold Standard (CPX/CPET – in this graph labeled as Ergometer).
- Seismofit® significantly overperforms vs. the substandard test methods.



# An experienced management team



Chief Executive Officer (CEO)

**Mikkel Kristiansen**

M.Sc. Investment/finance

- 10+ years experience from the MedTech industry, in various international leadership roles, focusing on commercialization and strategy



Chief Technology Officer (CTO), and co-founder

**Kasper Sørensen**

M.Sc. and PhD in Health Science and Technology

*Experience:*

- Research Assistant, AaU
- REDCap supporter at AaU



Chief Operations Officer (COO), and co-founder

**Peter Samuelsen**

M.Sc. Chem. Eng

*Experience:*

- 30+ years experience from the MedTech industry, including Coloplast as VP Corporate R&D and VP Corporate Venture
- Co-founder/ CEO/ Chair in Acarix A/S



Chief Medical Officer (CMO), and co-founder

**Peter Søgaard**

MD, Professor Cardiology

*Experience:*

- Associate Professor, Aalborg Uni. ,
- Co-inventor of multiple patents
- Co-inventor of SCG-fitness application



Chief Scientific Officer (CSO), and co-founder

**Samuel Schmidt**

PhD, Biomedical engineering

*Experience:*

- Associate professor, Aalborg Uni., and PhD supervisor
- Co-inventor of multiple patents
- Co-inventor of SCG-fitness application

# Some of our valued customers



# Some of our close partners



**AALBORG  
UNIVERSITET**



UNIVERSITY OF  
COPENHAGEN



Universität Hamburg  
DER FORSCHUNG | DER LEHRE | DER BILDUNG



**Sheffield  
Hallam  
University**

Advanced  
Wellbeing  
Research Centre



**Exercise  
is Medicine**



# Published clinical studies

## Validity and reliability of seismocardiography for the estimation of cardiorespiratory fitness

Mikkel T. Hansen, Tue Rømer, Amalie Højgaard, Karina Husted, Kasper Sørensen, Samuel E. Schmidt, Flemming Dela, Jørn W. Helge

Cardiovascular Digital Health Journal, 2023

DOI: <https://doi.org/10.1016/j.cvdhj.2023.08.020>

## Accuracy of a Clinical Applicable Method for Prediction of VO<sub>2</sub>max Using Seismocardiography

Mikkel Thunestvedt Hansen, Karina Louise Skov Husted, Mathilde Fogelstrom, Tue Romer, Samuel Emil Schmidt, Kasper Sørensen, Jørn Helge

International Journal of Sports Medicine. 2022

DOI: <https://doi.org/10.1055/a-2004-4669>

## A Chest-Mounted Accelerometer for Estimation of Cardiorespiratory Fitness

Samuel Emil Schmidt, Mikkel Thunestvedt Hansen, Kasper Sørensen, Tue Rømer, Lasse Gilemann, Dan Stieper Karbing, Mathias Krogh Poulsen, Jørn W. Helge, Peter Søgaard

2022 Cardiovascular Digital Health Journal – HRX San Diego

DOI: <https://doi.org/10.1016/j.cvdhj.2022.07.050>

## Estimation of Cardiorespiratory fitness using a using a chest mounted accelerometer

Samuel Emil Schmidt, Mikkel Thunestvedt Hansen, Tue Rømer, Peter Søgaard, Jørn W. Helge

2022 ESC congress – Barcelona

DOI: <https://doi.org/10.1093/eurheartj/ehac544.2769>

## Determination of Maximal Oxygen Uptake Using Seismocardiography at Rest

Mikkel Hansen, Birk Mygind, Tue Rømer, Kasper Sørensen, Samuel Emil Schmidt, Jørn W. Helge, Mathilde Fogelstrøm

2021 Computing in Cardiology (CinC)

DOI: <https://doi.org/10.23919/CinC53138.2021.9662756>

## Validity and reliability of a clinical non-exercise method for assessment of cardiorespiratory fitness using seismocardiography

M Thunestvedt Hansen, T Roemer, A Hoejgaard, K Husted, K Sørensen, S Schmidt, F Dela, J Wulff Helge

European Heart Journal, Volume 42, Issue Supplement\_1, October 2021, ehab724.3172

DOI: <https://doi.org/10.1093/eurheartj/ehab724.3172>

## A Clinical Method for Estimation of VO<sub>2</sub>max Using Seismocardiography

Kasper Sørensen, Mathias Krogh Poulsen, Dan Stieper Karbing, Peter Søgaard, Johannes Jan Struijk, Samuel Emil Schmidt

International Journal of Sports Medicine, 2020 sep.;41(10):661-668

DOI: <https://doi.org/10.1055/a-1144-3369>

# Further clinical studies on-going

- Royal Marsden Hospital (UK)
- Imperial College London (UK)
- Sheffield Hallam University (UK)
- St James's University Hospital, Leeds (UK)
- James Cook University Hospital, Middlesbrough (UK)
- York and Scarborough Teaching Hospitals (UK)
- Sheffield Teaching Hospitals (UK)
- Universität Hamburg (DE)
- German Sport University Cologne/Deutsche Sporthochschule Köln (DE)
- Southern Cross University (AU)
- Rigshospitalet (DK)
- Herlev Hospital (DK)
- Copenhagen University (DK)
- Aalborg University (DK)

# Seismofit® by VentriJect



simple  
fast  
reliable

estimation of  $VO_2$ max  
in less than 3 minutes

*“...probably the strongest indicator of health...”*



[www.ventriject.com](http://www.ventriject.com)

