

DETECTING SALMONELLA & BEYOND WITH N-LIGHT™ SALMONELLA RISK



SALMONELLA SPP. DETECTION – STATUS QUO

Salmonella infections pose a significant health risk, often linked to contaminated food products, making their control of food production crucial. The microbial environment in food production is extensive and intricate, with certain bacterial groups like Enterobacteriaceae and *Salmonella spp.* often at the centre of food safety concerns. Although *Salmonella spp.* is rarely found in factory settings, appearing in only 2 to 8 out of every 1000 routinely tested samples, total Enterobacteriaceae (TEB) counts can be much higher.

Unfortunately, no clear data links TEB counts to *Salmonella spp.* in environmental samples, which has been a significant challenge for researchers and industry experts. The EU Scientific Panel highlighted this issue and concluded that it is currently impossible to establish a direct correlation between Enterobacteriaceae and *Salmonella spp.*

! AYOPI GROUPS ARE AT THE HIGHEST RISK:



CDC estimates that *Salmonella* causes approximately 1.2 million illnesses and 450 deaths yearly in the US.



N-LIGHT™ SALMONELLA RISK – INDICATOR TEST

N-Light™ *Salmonella* Risk is not designed to detect only *Salmonella spp.* Instead, it also targets a small selection of closely related bacteria (such as *Citrobacter*, *Klebsiella* and *Enterobacter*) that share specific metabolic and genetic markers.

Positive N-Light™ *Salmonella* Risk does not mean finding exclusively *Salmonella spp.*


Instead, it highlights an elevated likelihood of contamination, allowing for timely interventions and corrective actions to ensure food production safety.

<i>Salmonella spp.</i>		Enterobacteriaceae
POSITIVE RESULT	POSITIVE RESULT	POSITIVE RESULT
<i>Salmonella spp.</i>	<i>Citrobacter</i> , <i>Enterobacter</i> , <i>Klebsiella</i> , <i>Salmonella spp.</i>	<i>Citrobacter</i> , <i>Enterobacter</i> , <i>Klebsiella</i> , <i>Salmonella spp.</i> , <i>Yersinia</i> , <i>Hafnia</i> , <i>Escherichia</i> , <i>Pantoea</i> , <i>Shigella</i> , <i>Serratia</i> , <i>Cronobacter</i> , <i>Samonia</i> , <i>Budvicia</i> , or <i>Proteus</i> , etc




WE TRANSFORM HOW PATHOGEN TESTING IS DONE

N-LIGHT™ SALMONELLA RISK TUBE CONTAINS A UNIQUE, INNOVATIVE TECHNOLOGY



Pioneering Phage Technology



Patented AquaSpark™ platform for Detection of Live Bacteria



The Biosafety Cap for Safe On-Site Enrichment

scan to see our technology in action




CERTIFIED AND TRUSTED

VALIDATED AGAINST
ISO 6579-1:2017


N-Light™ *Listeria monocytogenes*
Holds an AOAC® PTMSM
Certification



SIMPLE, EFFECTIVE TESTING PROCESS



SAMPLE



TRANSFER



INCUBATE





ACTIVATE





MEASURE


N-LIGHT™ SALMONELLA RISK TAILORED FOR MULTIPLE INDUSTRIES

- Fresh meat 
- Fermented dried meats

- Cheese/ raw milk 
- Milk powders
- Infant formula


- Nuts 
- Chocolate

- Spices 
- Egg powder

- Pet food 

A GROWING PORTFOLIO OF PATHOGEN DETECTION AND HYGIENE MONITORING

see other NEMIS products:



EMPOWERING YOU
TO FIGHT THE INVISIBLE
nemistech.com

