



Myth-busting gut pathogens

A microbiome-informed approach
to management



Introductions



Hannah Braye, MSc

Registered Nutritionist (BANT)
Clinical Director - Microba UK



A precision microbiome company delivering clinical-grade gut microbiome testing for healthcare professionals.



Founded in 2017 as a spin-out from the University of Queensland, Australia



Co-founded by Professors Philip Hugenholtz and Gene Tyson, renowned microbiome experts



Clinical-grade gut-health testing, using world-leading metagenomic sequencing technology




NATA-accredited laboratory, meeting ISO 15189 medical laboratory standards

First, do no harm.


HIPPOCRATES

Evolution in clinical thinking


Then:



Pathogen



Now:



Microbiome

Identify and eradicate pathogens and pathobionts with antimicrobials or antibiotics

Assess the whole microbiome and clinical case to address the underlying cause

Myth: All dysbiotic microbes are pathogens

Pathogen	a microbe that can cause disease e.g. <i>Entamoeba histolytica</i>
Opportunistic pathogen	a microbial strain that can cause disease in susceptible hosts e.g. <i>Klebsiella oxytoca</i>
Pathobiont	a microbial species associated with negative health outcomes e.g. <i>Bilophila wadsworthia</i>
Commensal	a microbial species associated with positive health outcomes e.g. <i>Faecalibacterium prausnitzii</i>



Myth: Because something natural, it is automatically gentler or more selective

- Action considered to be **selective** or **non-selective**.
- This refers to the specificity of the antimicrobial action.

Selective

Specifically suppress the growth of pathogens (e.g. *E. coli*)

Herbs: unknown due to limited research

Non-selective

May have antimicrobial effect on a wide variety of microbes including favourable commensal bacteria

Herbs: vast majority

Spotlight on berberine

An alkaloid found in several herbs and is often used to support blood sugar management

Study One by Zhang et al., 2020¹

- Multi-centre, double-blind, randomised, placebo-controlled study
- 409 newly diagnosed diabetic patients
- 7 days of antibiotics
- 600mg of berberine twice daily for 12 weeks
- Microbiome measured with metagenomics with FDR* corrections

1. Zhang et al. 2020, Nat Commun, doi:10.1038/s41467-020-18414-8

Study Two by Ming et al., 2021²

- Multi-centre, double-blind, randomised, parallel-controlled study
- 300 newly diagnosed hyperglycaemia patients
- 2-week run-in period with diabetes education
- 500mg of berberine twice daily for 16 weeks
- Microbiome measured with metagenomics without FDR corrections

2. Ming et al. 2021, Genome Med, doi:10.1186/s13073-021-00942-7

Increase in disease associated species

Bacteria	Disease Association and Features
<i>E. coli</i> ^{1,2}	Hexa-LPS producer High levels associated with Crohn's Disease and advanced liver fibrosis
<i>Klebsiella oxytoca</i> ¹	Hexa-LPS producer High levels associated with opportunistic infections and gastroenteritis
<i>Ruminococcus gnavus</i> ^{1,2}	High levels associated with IBS, Crohn's disease, atherosclerosis and obesity
<i>Klebsiella pneumoniae</i> ^{1,2}	Hexa-LPS producer High levels associated with hypertension, intestinal inflammation, NAFLD and Crohn's disease
<i>Bilophila wadsworthia</i> ¹	Hydrogen sulphide producer High levels associated with colon cancer and intestinal inflammation
<i>Klebsiella variicola</i> ¹	Hexa-LPS producer High levels associated with multiple myeloma and gestational diabetes
<i>Bacteroides_B dorei</i> ¹	High levels associated with colon cancer and type 1 diabetes
<i>Erysipelatoclostridium ramosum</i> ¹	High levels associated with obesity, type 2 diabetes, Crohn's disease, asthma
<i>Citrobacter koseri</i> ¹	Hexa-LPS producer
<i>Enterobacter cloacae</i> ¹	Hexa-LPS producer
<i>Enterobacter hormaechei</i> ¹	Hexa-LPS producer

Decrease in health associated species

Bacteria	Disease Association and Features
<i>Bifidobacterium longum</i> ¹	Reduction in harmful bacteria, anti-allergy effect, anti-obesity
<i>Bifidobacterium adolescentis</i> ¹	Anti-obesity effect, immune support Low level are seen in Crohn's disease, obesity, coeliac disease
<i>Bifidobacterium catenulatum</i> ¹	Low levels associated with colon cancer
<i>Bacteroides_F pectinophilus</i> ¹	Low levels associated with insulin resistance and metabolic syndrome
<i>Roseburia hominis</i> ¹	Primary producer of butyrate, supports immune system and reduces inflammation Low levels associated with IBD, hypertension
<i>Roseburia inulinivorans</i> ^{1,2}	Primary producer of butyrate Low levels associated with type 2 diabetes
<i>Roseburia intestinalis</i> ^{1,2}	Primary producer of butyrate Low levels associated with type 2 diabetes and Crohn's disease
<i>Ruminococcus_E bromii</i> ¹	Stimulates the growth of butyrate producing bacterial species
<i>Faecalibacterium prausnitzii</i> ¹	Butyrate producer
<i>Coprococcus eutactus</i> ¹	Butyrate producer
<i>Bifidobacterium longum</i> ¹	Reduction in harmful bacteria, anti-allergy effect, anti-obesity

Myth: You can tell if someone has a pathogen by symptoms alone

Acute pathogenic infection

- Watery diarrhoea, often sudden onset
- Urgency
- Blood or mucus in stools
- Abdominal pain and cramping
- Nausea and vomiting
- Fever
- Bloating
- Loss of appetite
- Fatigue / malaise

Chronic / unresolving infection

- Persistent or intermittent loose stools
- Steatorrhoea (greasy, floating stools)
- Unintentional weight loss
- Malabsorption & nutrient deficiencies (B12, iron)
- Persistent bloating and flatulence
- Chronic fatigue
- Post-infectious IBS-like symptoms
- Recurrent or relapsing episodes
- Reduced appetite


Test, don't guess

These symptoms overlap with many non-infectious conditions: IBS · IBD (Crohn's & ulcerative colitis) · coeliac disease · microscopic colitis · SIBO · bile-acid malabsorption · food intolerances · pancreatic insufficiency · medication effects · colorectal cancer

Myth: detection equals causation

The reality

- Many pathogens can be carried asymptotically, without causing disease
- Some sensitive detection methods can detect clinically insignificant levels of microbes
- The same finding can be clinically significant in one patient and irrelevant in another.
- Detection must be interpreted alongside symptoms and clinical context



Carriage rates vary widely with the pathogen, the detection method and the population sampled.

Finding a pathogen confirms it is present – not that it is causing the patient's symptoms

Myth: Every pathogen detected requires treatment

PATHOGEN MANAGEMENT GUIDELINES

Treat if detected



Entamoeba histolytica

Treat if symptomatic



Clostridium difficile pathogenic strains (toxin B; hypervirulent *C. difficile*); *Giardia lamblia*

Refer if symptomatic



Escherichia coli pathogenic strains (ETEC, EAEC, EPEC, O157, Shiga toxin, *Shigella* spp / EIEC)

Clinical judgement required — not necessarily pathogenic



Aeromonas, *Campylobacter*, *Cryptosporidium*, *Salmonella*, *Vibrio* spp.; *Cyclospora cayetanensis*; *Yersinia enterocolitica*; metagenomic potential pathogens (e.g. *Campylobacter_D upsaliensis*, *Clostridium_P perfringens*)

Pathogenic role unclear — exclude other causes of symptoms first



Dientamoeba fragilis; *Blastocystis* sub-types

Management depends on the specific microbe and the clinical presentation.

Indications for medical assessment

	Bacterial diarrhoea	IBD	Colo-rectal cancer
↑ Calprotectin	🚩	🚩	
↑ Lactoferrin	🚩	🚩	
Occult blood	🚩	🚩	🚩

Vulnerable patient

Immunocompromised

Elderly / Frail

Child

These indicate possible **severe or systemic infection**:

- Bloody diarrhoea
- High fever (>38.5°C)
- Severe abdominal pain or tenderness
- Persistent vomiting (unable to keep fluids down)
- Signs of dehydration
- Bruising easily or jaundice (suggests hepatic or haematologic involvement)
- Unexplained weight loss

Informed practice via whole ecosystem analysis



You need a complete picture of the microbiome and gut environment

DIVERSITY


What is the microbial diversity and richness?

SPECIES

Who makes up the whole ecosystem – bacteria, archaea, fungi, protists?

Which species are present and what is their abundance?

Which species are under or overabundant?



MICROBIAL FUNCTION

What is the whole microbiome's capacity to produce or consume metabolites associated with health?

GUT FUNCTION & ENVIRONMENT

Are there alterations in transit time, gut barrier integrity, digestive secretions or intestinal inflammation?

PATHOGENS

Are pathogenic bacteria or protist parasites present?

Microba combines multiple technologies for the most comprehensive assessment

Targeted PCR panel

Diagnostic pathogen
detection

Validated ELISA assays

Gastrointestinal health
markers

Metagenomic sequencing

Species

Diversity

Microbial function

Two complementary viewpoints

PCR



Targeted detection

Highly sensitive diagnostic detection of specific target pathogens

“Is this specific pathogen present?”

Metagenomics



Whole-ecosystem view

Analyses the whole microbial community — including pathogens you’re not specifically looking for, and their relative abundance

“What broader microbiome factors may be driving symptoms?”

+

Used together in the Microba test — a precise answer and the full picture.

Case Example - **Pathogen**

PCR Panel

Bacterial				
Detected Aeromonas spp. Bacterial	Detected Campylobacter spp. Bacterial	Not Detected Shiga Toxin Bacterial	Not Detected Enterotoxigenic E. coli... Bacterial	Not Detected Enteropathogenic E. coli... Bacterial
Not Detected Salmonella spp. Bacterial	Not Detected Enteroaggregative E. coli... Bacterial	Not Detected Yersinia enterocolitica Bacterial	Not Detected E. coli O157 Bacterial	Not Detected Vibrio spp. Bacterial
Not Detected Hypervirulent Clostridium... Bacterial	Not Detected Clostridium difficile toxin B Bacterial	Not Detected Shigella spp./EIEC Bacterial		
Parasitic				
Not Detected Entamoeba histolytica Parasitic	Not Detected Cryptosporidium spp. Parasitic	Not Detected Dientamoeba fragilis Parasitic	Not Detected Giardia lamblia Parasitic	Not Detected Cyclospora cayetanensis Parasitic

Metagenomic Analysis

Phylum	Species	Abundance	Prevalence	Distance from Average
 Campylobacterota	Campylobacter_D coli	0.40%	Rare	

Clinical context

Aeromonas spp.

- Food- and water-borne.
- Presents as asymptomatic carriage or traveller’s diarrhoea.
- Usually self-limiting
- Treat only if immunocompromised or symptoms are severe or persistent — guided by clinical presentation.

Campylobacter (C. jejuni & C. coli)

- Food-borne;
- Causes gastroenteritis.
- Usually self-limiting
- Treat only if immunocompromised or symptoms are severe or persistent — guided by clinical presentation.
- Faecal occult blood or suspected haemorrhagic colitis, arrange urgent investigation and specialist referral.

Gut Function & Environment

Gastrointestinal Health Markers

High Insights (8)

Calprotectin

Calprotectin is a protein released by neutrophils and is a non-specific marker...

Ref Range: < 50.00 µg/g



High Insights (4)

Lactoferrin

Lactoferrin is an iron-binding glycoprotein released from neutrophils and is a marker...

Ref Range: < 7.20 µg/g



Low

Faecal pH

Faecal pH is a marker of gut transit time. Low faecal pH may be indicative of rapid...

Ref Range: 5.70 – 7.30



Optimal

Secretory IgA

Secretory IgA (sIgA) is a marker of intestinal inflammation and increased intestinal...

Ref Range: 500.00 µg/g – 2000.00 µg/g



Optimal

Zonulin

Zonulin is a physiological regulator of intestinal permeability. Zonulin binds to a...

Ref Range: < 100.00 ng/mL



Optimal

Pancreatic Elastase

Pancreatic elastase is an enzyme produced by the pancreas that passes undigested...

Ref Range: ≥ 200.00 µg/mL






Not Detected

Faecal Occult Blood


Faecal occult blood is a marker of intestinal bleeding. Early diagnosis of faecal occult...

Diversity & Richness



Microbial Function

Microbiome Markers



Case Example - **Pathobiont**

PCR Panel

Bacterial

Not Detected

Aeromonas spp.
Bacterial

Not Detected

Enteropathogenic *E. coli*...

Not Detected

Hypervirulent *Clostridium*...

Not Detected

Yersinia enterocolitica
Bacterial

Not Detected

Clostridium difficile toxin B
Bacterial

Not Detected

E. coli O157
Bacterial

Not Detected

Enteropathogenic *E. coli*...

Not Detected

Campylobacter spp.
Bacterial

Not Detected

Salmonella spp.
Bacterial

Not Detected

Shiga Toxin
Bacterial

Not Detected

Vibrio spp.
Bacterial

Not Detected

Enterotoxigenic *E. coli*...

Not Detected

Shigella spp./EIEC
Bacterial

Parasitic

Not Detected

Giardia lamblia
Parasitic

Not Detected

Entamoeba histolytica
Parasitic

Not Detected

Cyclospora cayetanensis
Parasitic

Not Detected

Dientamoeba fragilis
Parasitic

Not Detected

Cryptosporidium spp.
Parasitic

Metagenomic Analysis

Q pathogen

×


All Associations

⌵

We've tested but no were detected


Gut Function & Environment

Gastrointestinal Health Markers



Microbial Function

Microbiome Markers



Most abundant species

All Species

Species	Phylum	Prevalence	Relative Abundance %	Distance from Average
Bacteroides_B dorei	Bacteroidota	Common	20.00 %	+2.62
Blautia_A wexlerae	Firmicutes_A	Very common	8.00 %	+1.57
Bacteroides uniformis	Bacteroidota	Very common	6.00 %	+1.81
Agathobacter rectale	Firmicutes_A	Common	6.00 %	+1.05
Lachnospira MIC8634	Firmicutes_A	Rare	6.00 %	
Faecalibacterium prausnitzii_K	Firmicutes_A	Common	3.00 %	+1.10
Ruminococcus_E bromii_B	Firmicutes_A	Common	3.00 %	+0.27
Blautia_A MIC9206	Firmicutes_A	Less common	2.00 %	+1.84
Lachnospira rogosae	Firmicutes_A	Common	2.00 %	+1.72
Faecalibacterium prausnitzii_G	Firmicutes_A	Very common	2.00 %	+0.94
Anaerostipes hadrus	Firmicutes_A	Very common	2.00 %	+0.73
Roseburia inulinivorans	Firmicutes_A	Common	1.00 %	+1.82
Bacteroides ovatus	Bacteroidota	Common	1.00 %	+1.71
Fusicatenibacter saccharivorans	Firmicutes_A	Very common	1.00 %	-0.97
Eubacterium_E hallii	Firmicutes_A	Common	1.00 %	+0.58
Faecalibacterium prausnitzii_C	Firmicutes_A	Common	1.00 %	+0.27
Blautia_A oheum	Firmicutes_A	Common	1.00 %	+0.89

Bacteroides_B dorei

This is a common inhabitant of the gut and is closely related to *Bacteroides vulgatus*.

Fuel Sources Used:

This species is a good degrader of fibre, a good degrader of mucin, and a moderate degrader of protein.

Metabolites produced:

Our genomic analysis indicates that most members of this species can produce the following metabolites: acetate, beta-glucuronidase, biotin (B7), branched chain amino acids, folate (B9), GABA, lactate, riboflavin (B2), vitamin K.

Metabolites consumed:

In addition, the genomic analysis shows that most members of this species can consume: GABA.

Emerging Research:

Higher levels of this species have been observed in patients with colon cancer and during the development of type 1 diabetes in children. This species has also been associated with diets high in red meat.

References

1

2





3

4

5

6

Richness and Diversity



Pathogen & Pathobiont Management Framework

Microba Microbiome Explorer Comprehensive

At Outset

Screen for **red-flag symptoms**. If present, medically investigate / refer.

Pathogen detected?

NO

Review other potential drivers of symptoms

e.g. Over/under abundant species, Functional dysbiosis, Gut function

YES

Non-resolving symptoms indicative of gut infection?

YES

ACTIVE PATHOGENIC INFECTION

Medical referral where relevant

Eliminate infection – Targeted antimicrobial therapy

Protect the microbiome – Prebiotics, *S. boulardii* / Bifidobacterium BB-12

Support the gut environment – Address post-infectious dysbiosis, intestinal permeability, intestinal inflammation, motility

NO

Consider risk factors:

microbe detected · patient vulnerability · red-flag GI markers

Investigation/ Intervention indicated?

YES

NO

PATHOBIONT OVERABUNDANCE

Avoid antimicrobial use

Competitive exclusion – e.g. GOS to target *E. coli*

Increase commensals – Dietary diversity, fibre, fermented foods, prebiotics (PHGG, RS2, beta-glucans)

Support the gut environment – gut barrier, motility, mucin production, digestive secretions

If diversity robust & overabundance severe / persistent despite supportive strategies:
Consider antimicrobials with simultaneous prebiotic actions (e.g. pomegranate husk)

Retest & Monitor

3–6 months

Key takeaways

1. Not all pathogens, opportunistic pathogens and pathobionts need eradicating – clinical context is key
2. Test don't guess - make informed clinical decisions by assessing whole microbiome health
3. Antimicrobials may not be the best intervention option for the gut microbiome and should be carefully considered
4. If anti-microbials are necessary, take steps to protect the wider microbiome



First, do no harm.

HIPPOCRATES

Designed to scale to clinical need

The tiered approach allows you to choose the right depth of analysis for each patient.

MICROBIOME EXPLORER COMPREHENSIVE

RRP £399

Metagenomics + GI Health Markers + Pathogen Screening

Full ecosystem mapping – A complete clinical picture from one sample.

MICROBIOME EXPLORER EXTENDED

RRP £349

Metagenomics + GI Health Markers

Microbiome and gut physiology insights – for chronic presentations once infection is ruled out.

MICROBIOME EXPLORER ESSENTIALS

RRP £299

Metagenomics

Focused microbiome analysis – Species, diversity and functional insights.



IPM Event Offer



50% off a Microbiome Explorer test

A comprehensive gut test with full microbiome ecosystem mapping and expert clinical interpretation.

USE CODE: IPM26_50

ORDER BY
30 JUNE 2026

microba.com Expect Better



REGISTER HERE

FIND US AT
BOOTH W37



Thank you for listening

Any questions?