ImmunoGuard PPE



Our Products

- 70% alcohol hand sanitiser
- Anti-bacterial wipes
- Sanitiser dispensers
- 3ply Type II R masks
- KN95 FFP2 masks
- Sterilised body gowns
- Sterilised body suits
- Anti-fog face shields
- Latex and nitrile gloves

Alcohol Hand Sanitiser

- Manufactured in Turkey
- 70% ethyl alcohol
- With alovera and vitamin E moisturising function
- Menthol fragrance
- PH 7.0
- Anti-bacterial, anti-viral and moisturising function
- Shelf life 3 years
- Bottle sizes 100ml, 150ml, 400ml, 1000ml, 5000ml
- ISO approved
- EU approved formula
- Shipment by road or air



Alcohol Hand Sanitiser





Alcohol Hand Sanitiser

Instructions for use: Dispense the required amount of solution onto dry palms to cover all areas of the hands. Rub the solution into the hands until the solution has dried from the hands. Do not use water or soap. Clean inside the nails and between fingers for best results. CAUTION: For external use only. Flammable liquid and vapour. Causes serious eye imitiation. In case medical advice is needed, bring the product container or label with you. Dispose contents / container in accordance with local / recional regulations.

Gebrauchsanweisung:

Zum hygienischen Händewaschen werden 2-3 ml des Produkts in die Handfläche gegeben. Mindestens 30 Sekunden. Bis die Hände und Handflächen verwendet werden. Wachen Sie die Innetsette der Nägel zwischen den Fingem und allen Innen-un Außenflächen.

Mode d'emploi :

Verser 3 mi sur mains propres et sèches. Frotter durant 30 secondes mains, poignets, espaces interdigitaux, pourtour des ongles. Ne pas nincer.

EN CAS DE CONTACT AVEC LES YEUX: Rincer avec précaution à l'eau pendant plusieurs minutes.

Précautions d'utilisation: A usage externe sur peau saine. Tenir hors de portée des enfants. Maintenir le récipient fermé de manière étanche. Tenir à l'écart de la challeur.

Modo d'uso: Versare 3mi sulle mani puite ed asciutte. Massaggiare 30 secondi, mani, polsi, spazi interdigitali, unghie. Non risciacquare. IN CASO DI CONTATTO CON GLI OCCHI: Sciacquare con acqua qualche minuto. Precauzioni d'uso Uso esterno: su pelle sana. Tenere fuori dalla portata dei bambini. Richuidere il flacone dopo l'uso. Tenere lontano da fonti di caiore.

INGREDIENTS FOR 100 ML	
Alcohol Denat %96	73,1 ml
Glycerin %85	1,45 ml
Aqua	25,05 ml
Aloe Barbadensis Leaf Extract	0,2 ml
Tocopheryl Acetate	0,2 ml









Avsedd användning:

Används för att rengöra händerna utan vatten och tvål.

Bruksarwisning:


Applicera gelen i handflatorna och grugga händerna mot varandra tills de känns torra och rena. Använd inte vatten eller tväl. FÖRSKTIGHET: Brandfarlig vätska och ånga. Orsakar allvarlig ögonirritation. Om medicinsk rådgivning behövs, ha produktbehållare eller etikett till hands. Förvaras oåtkomligt för barn. Kassera innehållet / behållaren i enlighet med lokala/regionala föreskrifter.

VARNING:

Fara/Lättantändlig/ Brandfarlig/Danger: Flammable

↑AKTUĞU KOZMETİK

EU Responsible Person: CPNP Reference: 3318309 CALIENTE GROUP LTD. Varna / BULGARIA 2: +90 541 317 63 13



Prd. : 01.04.2020 Exp. : 01.04.2023 Batch : L-122020







Anti-bacterial Wipes

- Manufactured in Turkey
- Anti-bacterial
- Kills 99.9% of germs
- Lemon fragrance
- PH 7.0
- Hydrogen peroxide sanitising effect
- Shelf life 3 years
- ISO approved
- CE marked
- Shipment by road or air





Hand Sanitiser Dispenser

- Manufactured in Turkey
- Rigid design
- High quality pump action
- Leak proof
- ISO approved
- Shipment by road or air








3ply Type II R Mask

- Fluid resistant
- CE, EN and FDA certified
- Filtration up to EN14683 standard
- BFE 99%
- ISO approved manufacturer
- Non-woven material
- 3 ply
- Dustproof, anti pollution, anti bacterial
- Sterile
- 175mm x 95mm
- Breathable and comfortable
- Shipment by airfreight





3ply Type II R Mask





KN95 FFP2 Mask

- Fluid resistant
- CE, EN and FDA certified
- Filtration up to EN14683 standard
- BFE 95%
- ISO approved manufacturer
- Non-woven material
- 3 ply
- Dustproof, anti pollution, anti bacterial
- Sterile
- 175mm x 95mm
- Breathable and comfortable
- Shipment by airfreight





Certification/ Class (Standard)	KN95 (GB2626-20 06)
Filter performance – (must be ≥ X% efficient)	≥ 95%
Test agent	NaCl
Flow rate	85 L/min
Total inward leakage (TIL)* – tested on human subjects each performing exercises	≤ 8% leakage (arithmetic mean)
Inhalation resistance – max pressure drop	≤ 350 Pa
Flow rate	85 L/min
Exhalation resistance - max pressure drop	≤ 250 Pa
Flow rate	85 L/min
Exhalation valve leakage requirement	Depressurizatio n to 0 Pa ≥ 20 sec
Force applied	-1180 Pa
CO ₂ clearance requirement	≤ 1%

KN95 FFP2 Mask







PPE Gown

- CE certified
- EN14126, EN13982, EN13034, EN1149 compliant
- Sterile
- Splash and flame proof
- ISO approved manufacturer
- Non-woven material
- One piece material, no stitching
- Knitted cuffs
- 50 gram
- 152cm x 80cm
- Breathable and comfortable
- Shipment by airfreight







PPE Gown















PPE Body Suit

- CE certified
- EN14126, EN13982, EN13034, EN1149 complian
- Sterile
- Splash and flame proof
- ISO approved manufacturer
- Non-woven material
- Laminated
- S, M, L, XL Sizes
- Breathable and comfortable
- Shipment by airfreight







PPE Body Suit

- Material 33gsm to 50gsm
- White
- Anti-static
- Sewing over locking and zipper
- Made from one piece
- Watertight zipper
- Chemical resistant





PPE Anti Fog Face Shield

- CE certified
- Anti-fog
- Splash and flame proof
- ISO approved manufacturer
- Reusable
- One size
- Breathable foam
- Comfortable
- Shipment by airfreight







Latex Gloves

- CE certified
- Powder free
- ISO approved manufacturer
- S, M, L sizes
- Comfortable
- Shipment by airfreight







Nitrile Gloves

- CE certified
- Powder free
- ISO approved manufacturer
- S, M, L sizes
- Comfortable
- Shipment by airfreight





