

# Drug & alcohol screening product range catalogue





At SureScreen, we believe that early, accurate diagnosis can significantly reduce the time to intervention, leading to better outcomes. Since our foundation in 1996, we have been providing proactive, bespoke solutions to organisations and healthcare providers.

We are an award-winning business that is known for high-quality products that meet and exceed our customer's expectations at a fair price.



## Contents

Detect	5
Drug Screening	8
Alcohol Screening	20
Confirmation Kits	26
Consumables	28
Contact Us	30
References	32

## Detect for your business

Protect your business, employees and customers with drug and alcohol screening to help ensure optimal workplace welfare.

Many industries are already faced with challenges to ensure employees' health and safety. Attributes like concentration, coordination and being of a sound mind when operating in a workplace environment are critical, not just to the success of a project, but for the safety of individuals and those working nearby. Incidents resulting from impairment or intoxication can be catastrophic and not only damage reputations but have the ability to stop a site, resulting in huge losses and potential slip on time-critical projects.

We can deliver a full end to end drug and alcohol screening program, as well as providing a comprehensive range of products and services to support you. We work with you to identify the best solution for your business:

- Review your policy and processes
- Identify the appropriate matrix, product and services
- Train your people to conduct screening to best practice
- Provide continual customer support

All our products and services are underpinned by our ISO9001 and ISO13485 accreditations, so you can be assured of the quality of our products and services.



## Windows of detection

The window of detection is the time that a drug or alcohol can be detected in a biological sample above a specified cut-off for the test being performed. The information provided is a guideline and does not indicate that a drug or alcohol will be detected for this long in all cases, as there are a number of factors that can influence the detection window.

### Urine

Substances can begin appearing in urine a couple of hours after a substance has been taken. These metabolites can be detectable in a urine sample for several days after use.

TYPE	DETECTION WINDOW
Drugs	12 hours - 12 days
Alcohol	Up to 24 hours

### Oral fluid

Oral fluid can be representative of the presence of a specified metabolite within the blood plasma. Therefore, it may indicate that the substance has recently been active in the body.

TYPE	DETECTION WINDOW
Drugs	Up to 72 hours

### Breath

Breath can be tested for alcohol using a hand-held device which provides immediate results.

TYPE	DETECTION WINDOW
Alcohol	Up to 24 hours



BZO	COC	THC
AMP	MEET	
OPI		


C T

SureScreen Diagnostics



Supra  
Lun PLY TPE  
Model: 500-02  
EN388: 2016  
EN12401: 2010  
EN12402: 2010






## Drug screening product range

Drug testing is the established way to quickly detect evidence of illegal drug use in a donor's sample. As drugs pass through the body they can be found in different matrices at different times.

Each drug test parameter looks for tiny particles of biological indicators of drugs using immunoassay technology. These particles are called metabolites, unique in their shape and form, they are highly distinguishable from each other. That's why our products provide reliable and accurate results.









## Rapid Cup

### Easy form of urine screening

The Rapid Cup is easy to use and provides test results available in 5 minutes. It is an efficient method of screening with a built-in temperature strip, privacy seal and multi-drug test.

### Benefits:

- ✔ Self-contained cup for collection and screening
- ✔ Individual strips for each drug tested
- ✔ Tight seal to prevent leaks
- ✔ Temperature strip to verify specimen validity
- ✔ Test for up to 18 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTS FOR
RCUP9DS1AB	9 Panel Rapid Cup	Urine	AMP/BZO/COC/K2/MDMA/MTD/OPI/PPX/THC
RCUP13DS6A	13 Panel Rapid Cup	Urine	AMP/BAR/BUP/BZO/COC/KET/MET/MTD/OPI/PCP/MDMA/TRA/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for urine testing products.



## Integrated Cup

### Controlled activation urine screening

The Integrated Cup is a multi-drug test and provides test results available at 5 minutes. It has a release key mechanism to control activation, along with a built-in temperature strip and privacy seal.

#### Benefits:

- ✔ Self-contained cup for collection and screening
- ✔ Tight seal to prevent leaks
- ✔ Temperature strip to verify specimen validity
- ✔ Activation by a key mechanism
- ✔ Test up to 12 parameters simultaneously
- ✔ Available with adulteration strip

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTS FOR
INCUP6DS1	6 Panel Integrated Key Cup	Urine	AMP/BZO/COC/MET/OPI/THC
INCUP6DS2	6 Panel Integrated Key Cup	Urine	AMP/BZO/COC/MTD/OPI/THC
INCUP6DS4300	6 Panel Integrated Key Cup	Urine	AMP/BZO/COC/MDMA/MOP/THC
INCUP7DS19300A	7 Panel Integrated Key Cup	Urine	AMP500/COC/MET500/MDMA/MOP/MTD/THC
INCUP10DS14	10 Panel Integrated Key Cup	Urine	AMP/BAR/BUP/BZO/COC/MET/MTD/OPI/PCP/THC
INCUP10DS13AB300	10 Panel Integrated Key Cup With Adulteration Strip	Urine	AMP/BAR/BUP/BZO/COC/EDDP/MET/MOP/PCP/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for urine testing products.



## Multi-Panel

### Ease of use in a clinical setting


The Multi-panel delivers ease of use for screening in a clinical setting, with test results available at 5 minutes. It utilises a multiline format to make efficient use of dip and read testing from a sample of urine across several parameters at the same time.

### Benefits:

- ✔ Cost effective
- ✔ Results can be photocopied for record-keeping
- ✔ Test up to 12 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTS FOR
S6DS1	6 Parameter Multi-panel	Urine	AMP/BZO/COC/MET/OPI/THC
S6DS2	6 Parameter Multi-panel	Urine	AMP/BZO/COC/MTD/OPI/THC
S6DS13	6 Parameter Multi-panel	Urine	AMP/BUP/BZO/COC/MTD/OPI
S6DS12	6 Parameter Multi-panel	Urine	BUP/BZO/COC/MTD/OPI/THC
S10DS2	10 Parameter Multi-panel	Urine	AMP/BAR/BZO/COC/MET/MDMA/MTD/OPI/TCA/THC
S10DS	10 Parameter Multi-panel	Urine	AMP/BAR/BZO/COC/MET/MTD/OPI/PCP/TCA/THC
S10DS1	10 Parameter Multi-panel	Urine	AMP/BAR/BZO/COC/MET/MDMA/MTD/OPI/PCP/THC
S12DS1	12 Parameter Multi-panel	Urine	AMP/BAR/BUP/BZO/COC/MET/MDMA/MOP/MTD/PCP/TCA/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for urine testing products.



## Multi-Cassette (Urine)

### Suitable for low-volume samples

The Multi-Cassette only requires a small amount of urine to perform the test, with test results available at 5 minutes. This multiline format test comes with a pipette for an easy-to-use process.

### Benefits:

- ✔ Multi-drug test incorporated
- ✔ Suitable for low-volume sample
- ✔ Test up to 12 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
S6MC1	6 Parameter Urine Multi-Cassette	Urine	AMP/BZO/COC/MET/OPI/THC
G10MC2300	10 Parameter Urine Multi-cassette	Urine	AMP/BAR/BZO/COC/MET/MDMA/MOP/MTD/TCA/THC
S12MC	12 Parameter Urine Multi-cassette	Urine	AMP/BAR/BZO/COC/MDMA/MET/MOP/MTD/OPI/PCP/TCA/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for urine testing products.



## Single Strip

### Designed for specific parameter testing

When monitoring or testing for one particular drug type, the simplest products to use are single test panels or single test strips, which are designed to look for the metabolite of a specific drug type.

### Benefits:

- ✓ Minimal cost solution
- ✓ Easy-to-use process
- ✓ Single parameter testing

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTS FOR
AMS	Amphetamine Single Strip	Urine	AMP
BAS	Barbiturate Single Strip	Urine	BAR
BZS	Benzodiazepine Single Strip	Urine	BZO
BUPS	Buprenorphine Single Strip	Urine	BUP
COS	Cocaine Single Strip	Urine	COC
CAS	Cannabis Single Strip	Urine	THC
MDMAS	Ecstasy Single Strip	Urine	MDMA
FENS	Fentanyl Single Strip	Urine	FYL
KETS	Ketamine Single Strip	Urine	KET
MTDS	Methadone Single Strip	Urine	MTD
MAMPS	Methamphetamine Single Strip	Urine	MET
OPS2000	Opiate Single Strip	Urine	OPI
SIDSPRE	Pregabalin Single Panel	Urine	PRE
TRAMS	Tramadol Single Strip	Urine	TRA

Note: Please refer to the references page for information on drug parameters and cut-off levels for oral testing products.



## SureSwab

### Ideal for workplace testing

The SureSwab is a rapid screening test to detect multiple drugs of abuse in oral fluid. It is the ideal solution for workplace or 'for cause' testing.

### Benefits:

- ✔ Large sponge to maximise sample collection
- ✔ Indicator to display sponge is fully saturated
- ✔ Strainer to extract sample and isolate collector sponge
- ✔ Tight seal to prevent leaks
- ✔ Concave base to deliver sample to test strips
- ✔ Hexagonal design can be photocopied for record-keeping
- ✔ Test up to 10 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETER TESTS FOR
SURESWAB6	6 Panel SureSwab	Oral fluid	AMP/BZO/COC/MET/OPI/THC
SURESWAB10	10 Panel SureSwab	Oral fluid	AMP/BAR/BUP/BZO/COC/MAMP/MTD/OPI/PCP/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for oral testing products.



## SureSwab Rapid

### Quick and easy form of oral drug screening

An easy-to-use screening test to detect multiple drug parameters in oral fluid. It is designed to be less invasive for improved donor experience and enables rapid collection time. A simple solution suitable for workplace, rehab and clinical testing.

#### Benefits:

- ✔ Small swab for easier saturation
- ✔ Indicator to display sponge is fully saturated
- ✔ Detection of both parent and metabolic Cannabis
- ✔ Slimline design for easy transportation or storage
- ✔ Test for up to 6 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
SURESWABRAPID6	SureSwab Rapid	Oral fluid	AMP/BZO/COC/MET/OPI/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for oral testing products.





## Twist Cassette

### Easy twist activation

The Twist Cassette is easy to use, simply twist the barrel to introduce the sample to the test strips, and read the results against the drug parameters printed on the test.

### Benefits:

- ✔ Multi-drug test incorporated
- ✔ Detection of both parent and metabolic cannabis
- ✔ Controlled activation by a twist mechanism
- ✔ Easy-to-use process
- ✔ Test up to 6 parameters simultaneously

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
ORAT6V3	6 Panel Twist Cassette	Oral fluid	AMP/BZO/COC/MET/MOP/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for oral testing products.


## Multi-Cassette (Oral)

Suitable for professional use

A Multi-Cassette including an individually wrapped collection sponge and collection vial that is ideal for professional bench-top testing. It uses multiline technology to test six drug parameters, making it appropriate for health professionals or laboratory use.

### Benefits:

- ✔ Minimal cost solution
- ✔ Easy-to-use process
- ✔ Multi-drug test incorporated




PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
ORAT6V1	6 Panel Oral Multi-Cassette	Oral fluid	AMP/COC/MET/OPI/MTD/THC

Note: Please refer to the references page for information on drug parameters and cut-off levels for oral testing products.



SureSwab

0  
Date  
Where



## Alcohol screening product range

Digital breathalysers are the gold standard for providing accurate and reliable results every time. With measurement in Breath Alcohol Content (BrAC) that goes beyond the legal drink drive limit, our product range helps monitor staff in safety-critical roles, understanding causation following an incident, or as part of a regular screening programme.

Our Alcohol screening products range from simple single-use breath screens to digital breathalysers.







## Alcometer

### Reliable and accurate alcohol meter

A breathalyser that can be used as an essential part of daily equipment for employers who have a zero-tolerance policy on alcohol.

### Benefits:

- ✔ Fuel Cell Sensor
- ✔ - 2.00mg/l BrAC detection range
- ✔ LCD 3-digit display
- ✔ Interchangeable mouthpieces for hygienic, multi-person testing
- ✔ Calibration period 6 months

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
ALCOM	Alcometer	Breath	ALC



## Intoximeter

### For increased durability

The Intoximeter offers legally defensible, quick and accurate screening for alcohol. As well as traditional breath testing, the product has a passive sampling cup and drink sniffer to detect alcohol in blown breath and outgassing of alcohol from suspect liquids.

### Benefits:

- ✔ Fuel Cell Sensor
- ✔ - 0.440g/210l BrAC detection range
- ✔ Backlit LCD display
- ✔ Interchangeable mouthpieces for hygienic, multi-person testing
- ✔ Calibration period of 12 months
- ✔ Passive sampling cup included
- ✔ Optional printer to provide evidential record
- ✔ Home Office approved

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETER TESTS FOR
INTOX	Intoximeter FST	Breath	ALC
INTOXWP	Intoximeter FST with printer	Breath	ALC



## Single Alcohol Test

### A cost-effective solution

For simple, rapid testing in breath, saliva or urine, these tests are perfect for on-the-spot testing. They can give a colour change that is read against a scale to determine the level of alcohol present.

### Benefits:

- ✔ Semi-quantitative results
- ✔ Cost-effective solution
- ✔ Easy-to-use process

PRODUCT CODE	PRODUCT NAME	MATRIX	PARAMETERS TESTING FOR
BSA	Breath Single Alcohol Test	Breath	ALC
SIDSETG1000	Urine Single Alcohol Test	Urine	ALC
AS2	Oral Fluid Single Alcohol Test	Oral fluid	ALC







## Confirmation Kits

Rapid screening allows you to eliminate the negative samples quickly and efficiently. If you do find a non-negative result, the next stage is to find out exactly what is in the sample and the quantity. This can be achieved by sending the sample to the laboratory where it is analysed using advanced chromatography methods.

Our confirmation kits make this process easy and include tamper-evident resources to ensure the legal validity of samples based on best practice chain of custody.



## Urine Drug Confirmation Kit

**Includes:**

- Consent form
- Tamper-evident bag
- Security seals
- Instruction sheet
- Absorbent pad
- Tamper-evident vials
- Pipette
- Disposable gloves

PRODUCT CODE	PRODUCT NAME	MATRIX
CONWC	Urine Drug Confirmation Kit	Oral fluid

## Oral Fluid Drug Confirmation Kit

**Includes:**

- Chain of custody form
- Tamper-evident bag
- Security seals
- Instruction sheet
- Saliva collectors
- Saliva swabs
- Disposable gloves


PRODUCT CODE	PRODUCT NAME	MATRIX
ORAT2CON	Oral Fluid Drug Confirmation Kit	Oral fluid

## Hair Drug Drug Confirmation Kit

**Includes:**

- Chain of custody form
- Tamper-evident bag
- Security seals
- Instruction sheet
- Hair specimen pouch
- Scissors
- Tweezers
- Piece of foil

PRODUCT CODE	PRODUCT NAME	MATRIX
HAK	Hair Drug Confirmation Kit	Hair



## Consumables



## Adulteration Strips

PRODUCT CODE	PRODUCT NAME	QUANTITY
S7ADLT25	Adulteration Strips	25 strips

## Mouthpieces for Breathalysers

PRODUCT CODE	PRODUCT NAME	QUANTITY
INTMP	Mouthpieces for Intoximeter	MOQ 100
ALCOMMP	Mouthpieces for Alcometer	MOQ 20

## Consent Form

PRODUCT CODE	PRODUCT NAME	QUANTITY
FORM.CONSENT	NCR Drug & Alcohol Collection Consent Form	25 strips

## Toilet Dyes

PRODUCT CODE	PRODUCT NAME	QUANTITY
DYE50/100	Toilet Water Dyes	Bag of 50 or 100
DYE	Individually Wrapped Toilet Water Dyes	1 piece

## Collection Cup

PRODUCT CODE	PRODUCT NAME	QUANTITY
CUP	Collection Cup	1 piece

## Disposable Gloves

PRODUCT CODE	PRODUCT NAME	QUANTITY
DG	Disposable Gloves	1 box of 50 pairs of gloves

## Contact us

Please don't hesitate to contact our Account Management team to find out more about SureScreen Diagnostics or our range of products. We can also discuss with you the bespoke products we are able to manufacture to suit specific requirements.

We are always happy to help, our specialist team is just a phone call or email away.

Call us: 01332 365 318

Email us: [info@surescreen.com](mailto:info@surescreen.com)

Visit our website: [www.surescreen.com](http://www.surescreen.com)



## References

The 'Reference Cut-Off Level' refers to the level that applies when the relevant drug code is referenced in a product description, unless otherwise stated. The table below provides information regarding the drug parameters available for urine testing.

The additional cut-off levels refer to alternative options available for bespoke product offerings. Please note parameters and cut-off levels vary by product.

### Drug Parameters and Cut-Off Levels for Oral Fluid Testing Products

FULL DRUG NAME	CODE	CUT OFF LEVEL (NG/ML)
Amphetamines	AMP	25/50
Barbiturates	BAR	50
Benzodiazepines	BZO	10/20
Buprenorphine	BUP	5
Cocaine	COC	15/20
Cotinine	COT	30/50
Fentanyl	FYL	10
Heroine	HRN	10
Ketamine	KET	50
Marijuana/Cannabis	THC	50
Methylenedioxymethamphetamine	MDMA	50
Methamphetamine	MAMP/MET	50
Morphine	MOP	40
Methadone	MTD	30
Opiates	OPI	30/40
Oxycodone	OXY	20
Phencyclidine	PCP	10
Propoxyphene	PPX	30/50
Synthetic Cannabis	K2	25
Tramadol	TML	30/50



### Drug Parameters and Cut-Off Levels for Urine Testing Products

FULL DRUG NAME	CODE	REFERENCE CUT OFF LEVEL NG/ML	ADDITIONAL CUT OFF LEVELS AVAILABLE
Amphetamines	AMP	1000	500/300
Barbiturates	BAR	300	200
Benzodiazepines	BZO	300	500/200/100
Buprenorphine	BUP	10	
Cannabis	THC	50	150/ 25
Clonazepam	CLO	300	
Cocaine	COC	300	100
Cotinine	COT	100	200
Ethyl Glucuronide	ETG	1000	
Fentanyl	FYL	1000	
Gabapentin	GAB	2000	
Hydromorphone	HMO	300	
Ketamine	KET	1000	500/ 300
(Lysergic Acid Diethylamide)	LSD	20	50
(6-Monoacetylmorphine)	6-MAM	10	
Methylenedioxymethamphetamine	MDMA	500	1000/ 300
Methamphetamine	MET/MAMP	1000	500/ 300
Morphine	MOP	300	100
Methylphenidate	MPD	1000	300
Methaqualone	MQL	300	
Methadone	MTD	300	200
Methcathinone	MCAT	1000	3000
Mephedrone	MEP	500	
Opiates	OPI	2000	
Oxycodone	OXY	100	300
Phencyclidine	PCP	25	
Pregabalin	PGB	2000	
Propoxyphene	PPX	300	
Synthetic Cannabis	K2	30	50
Tramadol	TML	100	
Tricyclic Antidepressants	TCA	1000	
2-Ethylidene-1, 5-Dimethyl-3, 3-Diphenylpyrrolidine	EDDP	200	300/100

Head office:

SureScreen Diagnostics,  
8B LittleOak Drive,  
Annesley, NG15 0DR

Tel: 01332 365318  
Email: [info@surescreen.com](mailto:info@surescreen.com)  
Web: [www.surescreen.com](http://www.surescreen.com)



SureScreen Diagnostics