

The RapidStop® Tourniquet



Available from
ProTrainings Europe Ltd
Phone 01206 809538
Click here for the online web

www.first-aid-online.co.uk



The RapidStop® Tourniquet relies on key technologies to achieve rapid, intuitive and easy one-hand self-application. Application is even faster when a responder uses both hands to apply the tourniquet to a victim.

The RapidStop® Tourniquet achieves full occlusion reliably and rapidly. It requires only gross motor control, allowing self-application to be readily possible. The RapidStop® Tourniquet is appropriate for controlling hemorrhage for injuries to both the arms and legs.

As Easy as 1-2-3 Three simple and fast steps set RapidStop® Tourniquet apart from the competition. Guided by instructions on the product, pocket-sized instruction cards and video format. Integrated instructions on the RapidStop® Tourniquet are included using stitched waterproof vinyl tags resistant to any wear and tear.

An individual who has sustained an injury to the arm or a leg and fully severed either the brachial, radial or ulnar artery has a window where stabilising the limb is essential to avoid the patient falling unconscious due to decreasing blood pressure. Achieving full occlusion in the shortest time possible is critical. Faster application means the injured are stabilised sooner and/or more victims can be saved.

1 PULL D-RING FIRMLY UNTIL SNUG



2 CRANK LEVER REPEATEDLY UNTIL BLEEDING STOPS



3 RECORD TIME, WRAP & TUCK EXCESS STRAP

