



CellAED™ TRAINER

LEARN. PRACTISE. REPEAT!

Designed to help you safely practise the application and use of your CellaAED®

Build confidence in using your CellaAED®

The first few minutes of a sudden cardiac arrest are critical. Knowing how to use your CellaAED confidently and acting quickly could be the difference between life and death¹.

That's why we've developed CellaAED *TRAINER*™ — a reusable replica of your CellaAED® so that you can safely practise over and over again.

Visit cormed-dx.com

1 St John Ambulance, Nov 2016

CorMed-DX



CellaAED® *TRAINER*™ App

Same size and feel as your CellaAED®

Visually differentiate your *TRAINER* device from CellaAED® by the yellow colour.*

Internal magnetic closure to repeat the snap activation

Push button to retract ribbon cable and reset the *TRAINER* to the start position

Reusable liner design covers the adhesive gel pads

Washable pads to practise applying in the correct position**

Comprehensive Training

Using your CellaAED *TRAINER*™ with the accompanying App is ideal for teaching and group training exercises to simulate real life scenarios of sudden cardiac arrest. The App enables you to hear the same audio voice prompts you would hear when activating and using your CellaAED® guiding you through the CPR and defibrillation process.

- Get familiar with hearing each voice prompt and practise the associated action at your own pace
- Play pre-set scenarios for infants and adults to gain confidence in positioning the pads correctly
- Choose from a variety of scenarios to test your listening and responsiveness and hone your skills in specific areas
- Access to instructional videos from identifying sudden cardiac arrest through to how to deactivate your CellaAED® after use

*For training purposes only and will not defibrillate

** Reusable pads can be removed and washed in warm soapy water to reinstate the adhesiveness as it wears over time

Visit cormed-dx.com

PROUDLY BROUGHT TO YOU BY



© 2023 CellaAED Life Saver Pty Ltd. All rights reserved. Unless otherwise noted all trademarks and the SNAP...PEEL...STICK trade mark are property of CellaAED Life Saver Pty Ltd.

DOC-3156v0