

Air-Q®3 Laryngeal Airways

**AVAILABLE
DIRECT FROM
FANNIN
OR
NHS SUPPLY
CHAIN**

Air-Q®3 Laryngeal Airways

The advanced, silicone Air-Q®3 helps deliver safety and confidence in airway management, no matter what airway challenges are faced. The unique ET tube ramp and epiglottis elevator enables swift and safe intubations.



Air-Q®sp3G

Uses positive pressure ventilation to self-pressurise the cuff. Cuff pressures cycle between peak airway pressure and PEEP levels. Features an extra wide gastric inlet that accommodate NG tubes up to 18Fr and two gastric channels.



Air-Q®3G

Features an extra wide gastric inlet that accommodate NG tubes up to 18Fr and two gastric channels.



Air-Q®sp3

Uses positive pressure ventilation to self-pressurise the cuff. Cuff pressures cycle between peak airway pressure and PEEP levels.




Air-Q®3

Provides exceptional ease of ETT placement when unexpected intubation is required.



Air-Q®3 Laryngeal Airway for perfect placement

- Anatomically correct placement of the masked airway
- Lifting of the epiglottis to maximise access
- Guide ramp to assist with ETT insertion and placement



Epiglottis elevator

The Air-Q®3 Laryngeal Airway includes an epiglottis elevator and raised heel to optimise airway placement and seal.





Air-Q®3 Laryngeal Airways – ordering details

Product code	Size	Ideal body weight (Kg)	Maximum OETT (mm)	Mouth opening (mm)	Tube length (cm)	Maximum gastric tube (Fr)	NHS Supply Chain code
Air-Q®sp3G							
60005	00	<2.0	3.0	5	6	5	FSM4368
60055	0.5	2-4	4.0	8	7	6	FSM4369
60105	1	4-7	4.5	11	9	8	FSM4476
60155	1.5	7-17	5	14	11	8	FSM4480
60205	2	17-30	5.5	17	14	10	FSM4475
60305	3	30-60	7	20	16	14	FSM4474
60405	4	60-80	8	23	18	16	FSM4473
60505	5	>80	9	25	20	18	FSM4479

For Air-Q®3G, Air-Q®3sp and Air-Q®3 product range codes, please contact your local Fannin representative