

Medical Technology | Creative Learning







CorMed-DX was established after 25 years of experience in providing NHS, Emergency Services, and Specialist Sectors with Technology and Services. With a Global network and a proven track record, we have the knowledge and connections to offer leading edge solutions and products to meet your needs.

Our team have a wide range of specialist skills:



medical technology



sales & marketing



tender bid/ management



learning



product development



large project roll out




TIME IS EVERYTHING™

Be prepared and confident to act quickly to save a life. CellAED® can be applied in seconds.

Graphics show how to apply the device on both adults and infants

A perforated label marks the place for easy snap activation Audio speaker gives instructions on what to do every step of the way



Infant mode button for use on infants weighing less than 10kg Length: 19.6cm Width: 9.3cm Depth: 1.7cm

Weight: Approx 300g

LED indicator with periodic blinking showing the status of the device

Your handheld, smart, personal defibrillator

CorMed-DX have partnered with Rapid Response Revival as the Authorised Resellers of CellAED® in the UK.

CellAED® is designed to save lives from sudden cardiac arrest, which can happen to anyone, anywhere, at any time.

During a cardiac arrest, every minute without an AED lowers the chance of survival by 10%!

With CellAED®, you can ensure there is always an AED nearby to help save a life.

Visit: www.cormed-dx.com



Power CellAED® by snapping in half along perforated centre line. This will automatically start the audio instructions.



Remove the protective liner under your CellAED®. This will expose the integrated gel pads and extend the connecting cable.

PEEL



For adults, place each half of the CellAED® on the bare chest. For infants, place one half on the infant's chest and the other half on the back.

CorMed-DX




tac:life c

dynamic.



LEARN. PRACTISE. REPEAT.



DESIGNED TO HELP YOU SAFELY PRACTICE THE APPLICATION AND USE OF YOUR CellAED®

Practice using your CellAED TRAINER™ with the same audio voice prompts included in your CellAED®, which guide you through the CPR and defibrillation process.

- Get familiar with hearing each voice prompt and practise the associated action at your own pace
- Play pre-set scenarios for infants and adults to gain confidence in positioning the pads correctly
- Choose from a variety of scenarios to test your listening and responsiveness and hone your skills in specific areas
- Access to instructional videos from identifying sudden cardiac arrest through to how to deactivate your CellAED® after use

Visit: www.cormed-dx.com

CellAED® Trainer™ includes a free APP which gives scenarios and audio prompts so you and your teams can practise how to use the device in an emergency on adults or infants.



CorMed-DX







tac:life

respond first with confidence



Covert and overt options for military, medical, close protection, and civilian use.

Thanks to the IFAB™ and CellAED®, it is now possible to carry an AED, as well as life-saving bleed control kit including the 'Big Four' you'll need to control traumatic haemorrhaging wherever you go.

The first-of-its-kind Individual First Aid Belt (IFAB™) gives users unparalleled peace of mind knowing that they have the core equipment required to save a life on their person, at all times.

The IFAB™ is made up of two main components:

THE MAGCORE™

the mounting platform for the IFAB. It has two metal plates that line up with reciprocating magnets on the second component, the MAGCARD™. It is low profile and sits comfortably on the user's lower back/waist.

The MAGCORE™ comes in two variants:

- The MAGCORE Molle™ attaches to any existing molle pals retention system.
- . The MAGCORE Belt™ attaches to any belt up to two inches in width.

THE MAGCARD™

attaches to the MAGCORE via its bespoke rare earth magnets with more than 12kg pull, and velcro. Its pouches and sleeves, made of the highest quality materials, carry the AED and all the kit you need to control the bleed.

The IFAB AED™

with full control bleed kit weighs less than a kilogram (dependant on products included), and measures 31cm by 12.7cm.

Included in the IFAB AED™: Tourniquet, Chest Seal, Basic Life Support Card, Trauma Shears, Haemostatic Dressing, 4" Compression Bandage, Marker Pen, and a CellAED (Razor - optional extra).



CorMed-DX



tac:life



F+RSTAED

Responder





The FirstAED responder box is designed to provide all the kit and equipment needed to tend to trauma and cardiac incidents. Complete with a bleed control kit and the revolutionary CellAED®, all in one place, fit for purpose.

The box can be mounted on walls as well as transported in vehicles or carried. This equipment is an affordable and versatile option to make sure your family remains safe. This can be wall mounted, carried, or stored in a vehicle, making sure you have the right kit there when you need it.



Visit: www.cormed-dx.com

The FirstAED Responder Box Includes...

- Items held in waterjet cut foam inserts
- 2 x CAT tourniquets
- 2 x Trauma Bandage 4"
- 2 x Foxseal chest seals
- 2 x Celox rapid
- 2 x pair of nitrile gloves
- 1 x CellAED®
- 1 x Tuff cuts
- 1 x Gallant disposable razor blade
- 2 x NPA (6 & 7)
- 1 x Sharpie pen
- 1 x ATMIST/MARCH card

CorMed-DX



tac:life

dynamic.

dynamic.

creative learning



As one of the UK's leading elearning design companies, dynamic. creates and delivers high-impact online learning for your audience. It's a highly effective way to drive efficiency, knowledge, productivity profitability throughout your organisation. Whether you need bespoke content, off-the-shelf elearning or you'd benefit from our unique dynamicLMS, you can discuss your options and opportunities on how you can add value to your business.

solutions

Using the latest creative and immersive technology, dynamic. offers a wide range of solutions, including:



Advanced Authoring




Games



360° Virtual Tours



VR (Virtual Reality)



AR (Augmented Reality)



Apps & other support materials

CorMed-DX



Emergency Services

the dynamic approach

Dynamic have extensive experience of creating comprehensive materials to support those working in both health and social care sectors. We are keen to look at where our expertise could be beneficial to those who work within Emergency Services and have taken time to consider some of your potential hotspots:



Medical Training

dynamic. have created materials to support many medical procedures. Whether you have a need for ad-hoc bespoke courses or programmes to support your mandatory learning agenda dynamic. can help.

Vehicle Layout

We appreciate that the effective configuration of an ambulance is essential. dynamic. could support with the development of materials to help with training on the clean-down and make-ready processes and vehicle safety checklists.

Bespoke Learning



We have the skills to create immersive materials that draw upon the use of VR for training. The example shown is based on fire safety training provided for those working in acute hospital setting.



We can develop materials that put the learner at the centre of a scenario. The example shown relate to resuscitation in neonatal care and scenarios in a socail care setting.



The Scan and Learn technology helps to put learning at the point of need. You can have cards or even products that can be scanned, realtime, to display instruction or guidance. This can include videos, tutorials, and even links to purchase a product. The example shown is the CellAED® Scan and Learn App, simply point the camera at the CellAED® or CellAED Trainer® to reveal the guidance: Instructional Video, CellAED Trainer® App download, CellAED® online shop link.

Scan and Learn

CorMed-DX



tac:life dynamic.







Visit cormed-dx.com

Email:

info@cormed-dx.com

Phone:

0800 3317999

Visit: www.cormed-dx.com

Address:

CorMed-DX Ltd.

Aizlewood's Mill,

Nursery Street,

Sheffield,

S3 8GG

CorMed-DX





tac:life dynamic.