

PRESS RELEASE

Accora Announces Release of FloorBed for Plus Size Service Users

When a single profiling bed is too narrow and a bariatric profiling bed is just too wide....

CHALLENGE:

Single bed too narrow

inadequate turning space for wider clients
inadequate space for those with contractures
inadequate space for those used to sleeping in a double bed

Bariatric bed too wide

inadequate space in environment
inability to reach bed levers
increased risk to carers having to reach over the width of the bed to manage care tasks.

SOLUTION:

Accora are proud to announce the launch of FloorBed1 Plus. At 105cm wide, FloorBed1 Plus is 15cm wider than a standard single; the extra width in conjunction with the ultra-low function of 7cm the FloorBed Plus will meet a wide variety of client's needs. In addition, the FloorBed Plus will:

- ✓ Reduce a bed fall to a roll with the aim to minimise injury
 - ✓ Enable independent leg lifting
 - ✓ Promote independent
- ✓ bed transfers and reduce carer requirements

Height range: 7 – 64cm

Width: 105cm

Weight capacity: 35stone

- ✓ The additional 15cm width gives extra space which may be required when turning
- ✓ The extra width is particularly beneficial when there is not adequate space for the full width of a bariatric bed or if carers would struggle to reach across a bariatric bed to carry out care tasks
- ✓ At 105cm wide this bed may reduce the risk of bed falls for a person who is used to sleeping in a single bed

Accora will be showcasing this product at the [OT Show 2016](#), Stand K40. To find out more about appropriate applications please book a seat at one of the [sessions](#):

23rd November 2016 @ 11am or 24th November 2016 @ 2pm

Take a sneak preview with our short video: [FloorBed1Plus](#)