

X Series® Advanced

ZOLL®



MOVING CARE FORWARD

THE GROUNDBREAKING NEW MONITOR/DEFIBRILLATOR FR





OM ZOLL

To improve patient outcomes, healthcare providers rely on feedback and data at every step of the rescue. Access to real-time clinical feedback on the scene can take your care to the next level and enable you to manage multiple disease states.

Real-Time Clinical Feedback—Available Only Through ZOLL

For years, ZOLL® devices have provided real-time clinical feedback on CPR to help rescuers deliver high-quality compressions to patients in cardiac arrest. The X Series® Advanced monitor/defibrillator goes beyond cardiac arrest, providing real-time clinical feedback for patients in respiratory distress and patients with a traumatic brain injury:



Real CPR Help® – provides real-time clinical feedback on compression rate and depth



Real BVM Help™ – provides real-time clinical feedback on ventilation volume and rate



TBI Dashboard™ – provides the right information at the right time to help rescuers effectively manage traumatic brain injury (TBI) patients


DISEASE STATE MANAGEMENT AND ACCESS TO KEY DATA

Manage Multiple Disease States

In the pre-hospital environment, access to data helps you manage trauma, cardiac arrest, chest pain, respiratory distress, pediatric issues, and other conditions. X Series Advanced provides information and feedback on CPR, ventilation, and TBI, increasing the range of treatable patients when compared to other currently available monitor/defibrillators.




Standard Monitor/defibrillator



- Cardiac Arrest
- STEMI
- Other

X Series Advanced



- Cardiac Arrest
- STEMI
- Trauma
- Respiratory
- Other



Seamless Device Integration

For cardiac arrest patients, the X Series Advanced features Real CPR Help and Real BVM Help, providing simultaneous real-time feedback on both CPR and ventilation quality. When used with the AutoPulse®, ZOLL's automated CPR device, the X Series Advanced captures the unique chest compression data and displays it on screen.

Easy Access to Data

With RescueNet® CaseReview, clinicians can view data from the X Series Advanced post-case for effective QA/QI—no more chasing crews to get case files. With the push of a button, medics can send files directly to your secure system for performance-focused debriefing. Access to data, including ventilation and CPR performance, is as fast as opening a web browser.






NEED MORE
INFORMATION?

Please contact your local ZOLL representative or distributor/partner, or email ServiceContractsAdmin@zoll.com

The Power of Telehealth

With Remote View, ZOLL's telehealth technology, healthcare professionals can view data from the X Series Advanced from a remote location. Remote View provides effective decision support that enables EMS providers to meet their telehealth needs.

X Series® Advanced Transmission



Cellular Tower



Remote View Server



Web Browser





ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | contact your local or partner | www.zoll.com

Copyright © 2020 ZOLL Medical Corporation. All rights reserved. X Series, and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners. ZOLL X Series Advanced is not available for sale in the United States. The product has not received regulatory clearance/approval by the Food and Drug Administration.
Printed in U.S.A. MCN EP 2006 0314-05

ZOLL®