

TITAN TI TIROL STRETCHER & KIT

Part of the stretcher range from Lyon Equipment


The TITAN TI SPLIT TAPERED STRETCHER and TIROL KIT is a lightweight and versatile patient transport system ideal for long haul retrieval and technical rescue on even the most difficult terrain.

The TIROL KIT is designed to fit a purposebuilt TITAN TITANIUM SPLIT BASKET STRETCHER which has additional head and foot plates incorporated. The combination of these two lightweight products makes it easier to operate in hostile and inaccessible areas.



See overleaf for detailed product features

CODE	DESCRIPTION	WEIGHT	CONFORMITY
II-0II9TIROL	TITAN TI SPLIT TAPERED STRETCHER with handle mounting plates	9 kg	UKCA, CE, Class Medical Device
TIROL-KIT	TIROL KIT includes light Mule, wheel system, handles, shoulder supports, storage bag	10.8 kg	N/A





At Lyon Equipment, we design, manufacture and supply high quality and innovative, adventure, work at height and specialist rescue equipment - delivering expert advice on all our products and their applications



TIROL features:

- Split apart design for ease of storage and transportation
- Lightweight Titanium frame
- Built in head and foot attachment points for handles
- 25.4 mm top tube allowing rescuers to maintain a good grip, especially in long carries
- Four StratLoad[™] attachment points that make it quick to clip your stretcher bridle
- UHMW-PE Back Support
- High Density Polyethylene netting -HDPE

TIROL KIT features:

- The kit is lightweight and can be split between rescuers
- Adjustable handles to suit the rescuers and terrain
- Handles can be positioned to function as legs for the structure and provide a stationary position for the stretcher
- Lightweight single wheel system including fast-clip locking, designed for mountain terrain
- 225 kg capacity with SWR safety factor of 1.5 to guarantee secure performance of the product

TIROL KIT includes:

- Lightweight single wheel system
- 4 Adjustable carry handles
- 2 Shoulder supports
- Storage bag

Further resources:

Contact us for further information on servicing and repairs

Visit www.lyon.co.uk/videos for short how to guides on how to-

- Split the stretcher into 2 sections
- Assemble the TIROL KIT WHEEL
- Use the TIROL KIT HANDLES

Product accessories:



Lyon Equipment offer other stretcher solutions -

To learn more about our products and the services we offer please visit www.lyon.co.uk

Lyon Equipment LimitedRegistered Office - Unit 3-7, Tebay Business Park, Old Tebay, Penrith, Cumbria, CA10 3SS, United Kingdom Company number 3107118 - Registered in England & Wales.

015396 24040

WRsales@lyon.co.uk

www.lyon.co.uk

