




Is your defib emergency ready?



In a cardiac arrest, every second is vital.

Early CPR and defibrillation can more than double the chances of surviving an out-of-hospital cardiac arrest, but if defibs aren't visible to ambulance services, how can they direct people to them in an emergency?

Registering your defib on The Circuit – the national defibrillator network, ensures ambulance services can locate your defib when 999 is called.

It is free and takes as little as 5 minutes.



Find out if your
defib is registered

Over 100,000 defibs are registered and ready to save a life, but is yours?

Find out at defibfinder.uk

**Not registered? Register
today at thecircuit.uk**

The Circuit helps ambulance services direct bystanders to their nearest defib, so defibrillation can start as early as possible.

Are you the guardian of a defib that is yet to be registered on The Circuit?

Once registered, you will receive reminders whenever something needs to be done – such as status checks or replacing expired pads. By regularly updating your defib's status on The Circuit, ambulance services can trust it's ready to use in an emergency.

Next steps

- 1 Visit **thecircuit.uk**
- 2 Set up your account
- 3 Register your defib
- 4 Regularly update your defib's status on The Circuit



Scan the QR code to go straight to registration.

Information to have to hand

The location and surroundings of your defib

The brand, model and device serial number

When the defib is accessible

When the pads expire and if it is emergency ready

Brought to you by




In partnership with



Resuscitation Council UK



Proudly supported by
 omaze

In 2025 Omaze raised funds for British Heart Foundation that will help to power The Circuit for two years.

TheCircuitLft2025