



ProVu™ Mini Video Laryngoscopes

Introducing our new
Mini Handle Paediatric Blades
to the ProVu Video
Laryngoscope range.


Slimline shape yields a more tactile feel
for the user, ideal for when dealing
with delicate anatomies.

Combines direct laryngoscope techniques
with video technology.



100% metal blades.
Provided sterile.
Single Use.

Standard handles and
Adult blade sizes also available.




Miller Blades

An improved view with traditional
laryngoscopy techniques.



Macintosh Blades

Low profile blade, ideal for
smaller mouth openings.



Hyperangulated Blades

Difficult airway blade for
improved laryngeal exposure.



Designed with the most vulnerable patients in mind

Use of oversized handles may compromise intubation success, and carries a greater risk for injuring the oropharynx.¹

Slimline shape, lightweight and ergonomically balanced makes ProVu Mini Handle ideal for use with smaller patients.


Ultimate Visualisation

Get everyone involved.
ProVu mini Blades can be connected to the ProVu 8" display, for improved visibility to you and the team.

8" Display + Pole Stand and Articulating Arm

Achieve the optimal screen viewing angle for you and the team when mounted onto the ProVu articulating arm.

Storage for ProVu devices in the basket and charging capabilities on the pole.



Single Use ProVu Video Laryngoscope Blades

STANDARD HANDLE	MINI HANDLE	DESCRIPTION	QTY
040-08-0110	040-09-0110	ProVu Single Use Video Laryngoscope Handle with Macintosh Blade Size 1	10
040-08-0120	040-09-0120	ProVu Single Use Video Laryngoscope Handle with Macintosh Blade Size 2	10
040-08-2110	040-09-2110	ProVu Single Use Video Laryngoscope Handle with Hyper Angulated Blade Size 1	10
040-08-2120	040-09-2120	ProVu Single Use Video Laryngoscope Handle with Hyper Angulated Blade Size 2	10
	040-09-0410	ProVu Single Use Video Laryngoscope Handle with Miller Blade Size 00	10
	040-09-0400	ProVu Single Use Video Laryngoscope Handle with Miller Blade Size 0	10
	040-09-0401	ProVu Single Use Video Laryngoscope Handle with Miller Blade Size 1	10

Reusable Displays and Mounting Accessories

Adult sizes also available.

PART	DESCRIPTION	QTY
040-07-0035	ProVu 3.5" Display + Charging Dock	1
040-07-0080	ProVu 8.0" Display + Charging Cable + 2m Amplifier Cable	1
038-990-101	ProVu Display Pole Stand + Universal Articulating Arm	1
038-990-104	ProVu 8.0" Display Universal Articulating Arm	1
038-990-102	ProVu Display Table Stand	1

Reusable Chargers and Cables

040-07-0020	ProVu Universal Reusable 2m Amplifier Extension Cable	1
038-990-013	ProVu Dual Micro USB 1.0m Charging Cable	1
038-990-010	ProVu Charging Dock for 3.5" Display	1

¹Christine E Whitten, MD - 10 Common Pediatric Airway Problems- and their solutions - Anesthesiology news airway management 2019.

Flexicare Medical Limited

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK
+44 (0)1443 474617 • enquiries@flexicare.com • www.flexicare.com

FC7682C | 122023

